

John P. Wilson^{1*}, Visa Meyyappan², Domenic N. Narducci², Ben A. Neely³, Jim A. Laugharn², Darryl J. Pappin^{1,4}

IFITM

*To whom correspondence should be addressed: john@protifi.com



Covaris®

Component	Variance (%)
Extraction	72
Instrument variance	16
Instrument stability	8
Digestion	3

SDS/AFA-S-Trap workflow for FFPE samples.

Fig 8 (left): Average coefficient of variation (CV) of protein extraction yields across all experiments. Coefficients of variance (CVs) of AFA/SDS extraction were consistently $< 10\%$ ($5\% \pm 3\%$).

Organ	Condition	Cytoplasm	Membrane	Nucleus
Brain	1% 500Iop-iod WtEA AHA	1850	1750	1150
	12% 125-iod-WaterCSD AHA	1150	850	750
	10% 125-iod-WaterCSD(25% acetone) AHA	900	700	600
	TP AHA	1800	1650	1050
	10% 125-iod-WaterCSD(25% acetone) SB	850	900	750
Pancreas	10% 500Iop-iod WtEA AHA	1450	1350	1050
	12% 125-iod-WaterCSD AHA	850	850	550
	10% 125-iod-WaterCSD(25% acetone) AHA	650	650	550
	TP AHA	1450	1300	1050
	12% 125-iod-WaterCSD(25% acetone) SB	650	450	450

The top diagram shows frequency in Hz (1 Hz, 1 kHz, 1 MHz) on the x-axis, with regions for infrasonic, sonic, and ultrasonic. It highlights 'Sonic radar' in the sonic region and 'Diagnostic imaging' and 'Covaris' in the ultrasonic region.

The bottom diagram shows wavelength in nm (1m, 100nm, 10nm, 1mm, 100um, 10um) on the x-axis, with regions for 'Sonic' and 'Ultrasonic'. It highlights 'Sonic radar' in the Sonic region and 'Diagnostic imaging' in the Ultrasonic region.

The diagram illustrates the S-Trap workflow in five steps, connected by arrows and a central time indicator:


- Step 1:** A microcentrifuge tube containing a grey liquid. Below it, the text reads: "Protein in 5% SDS (after AFA or)".
- Step 2:** An arrow points to a tube with a green bead layer and a grey liquid layer. Above the arrow is the text: "Add S-Trap™ binding buffer, load".
- Step 3:** An arrow points to a tube with a green bead layer and a clear liquid layer. Above the arrow is the text: "Spin, wash, spin".
- Step 4:** An arrow points to a tube with a pink bead layer and a clear liquid layer. Above the arrow is the text: "Add protease e.g. trypsin". Below the tube is the text: "Protein trapped, denatured, cleaned".
- Step 5:** An arrow points to a tube with a pink bead layer and a pink liquid layer. Above the arrow is the text: "Digest (1 hr 47 °C)".
- Step 6:** An arrow points to a tube with a pink bead layer and a pink liquid layer. Above the arrow is the text: "Spin elution". Below the tube is the text: "Digested peptides".

A large double-headed arrow at the top spans from the second step to the sixth step, with the text "A few minutes" above it.


96-well plate

Midi Mini Micro

Condition	% known brain-specific proteins observed
500/50 mM TEAB AFA	51%
125 mM NH ₄ HCO ₃ AFA	31%
125 mM NH ₄ HCO ₃ /20% TP AFA	26%
TP AFA	50%
125 mM NH ₄ HCO ₃ BB	30%




	Extraction buffer	Extraction technique	Digestion technique	Total IDED proteins	% of proteins IDED in any condition
Blue	SDS	AFA	S-Trap	3887	82.5%
Red	NH ₄ HCO ₃	AFA	S-Trap	3779	80.2%
Green	SDS	BB	S-Trap	3276	67.4%
Yellow	NH ₄ HCO ₃	BB	S-Trap	2211	46.9%



	Extraction buffer	Extraction technique	Digestion technique	Total IDED proteins	% of proteins IDED in any condition
Blue	SDS	AFA	S-Trap	3091	83.3%
Red	NH ₄ HCO ₃	AFA	S-Trap	2307	62.2%
Green	SDS	BB	S-Trap	2234	60.2%
Yellow	NH ₄ HCO ₃	BB	S-Trap	1813	48.9%

	Extraction buffer	Extraction technique	Digestion technique	Total IDED proteins	% of proteins IDED in any condition
Blue	SDS	AFA	S-Trap	2508	97.0%
Red	NH ₄ HCO ₃	AFA	S-Trap	52	2.0%
Green	SDS	BB	S-Trap	581	22.5%
Yellow	NH ₄ HCO ₃	BB	S-Trap	94	3.6%



1. Sources of Technical Variability in Quantitative LC-MS Proteomics. Piewhowsky et al. *J. Proteome Res.* 12(5) 2013.
2. Laugharn JR, J.A. and Durin, G., Covaris Inc, 2016. Method and apparatus for shearing of genomic material using acoustic processing. U.S. Patent 9,486,756. See also <https://covaris.com/resources/publications/>.
3. Zougman, A., Selby, P.J. and Banks, R.E., 2014. Suspension trapping (S-Trap) sample preparation method for bottom-up proteomics analysis. *Proteomics*, 14(9), pp.1006-1000.
4. Ludwig, Katelyn R., Monica M. Schroll, and Amanda B. Hummon. Comparison of In-Solution, FASP, and S-Trap Based Digestion Methods for Bottom-Up Proteomic Studies. *Journal of proteome research* (2018).
5. HaileMariam, Milkesha, Rodrigo Vargas Eguetz, Harinder Singh, Shiferaw Bekele, Gobena Ameni, Rembert Pieper, and Yanbao Yu. S-Trap is an ultrafast sample preparation approach for shotgun proteomics. *Journal of proteome research* (2018).