

Recommendations for Different Proteases

S-Trap's™ work with your choice of protease.

Protease	Protease:Protein (w/w)	Digestion Time	Digestion Temperature	Digestion Buffer	Digestion Buffer for Isobaric Labeling
Asp-N ¹	1:20	3 hours	37 °C	0.5 mM ZnCl ₂ in 50 mM NH ₄ HCO ₃	0.5 mM ZnCl ₂ in 50 mM TEAB
Chymotrypsin ²	1:30	16-18 hours	37 °C	50 mM TEAB, 1 mM CaCl ₂	
Glu-C ³	1:10	16 hours	37 °C	0.1 M Phosphate Buffer, pH 7.7*	
Lys-C ⁴	~1:30	4 hours	37 °C	Tris-HCl Buffer, pH 8.5	50 mM TEAB
Trypsin	1:10	16-18 hours	37 °C	50 mM TEAB, pH 8.5	

*Sodium phosphate must be removed prior to MS-analysis, especially when performing phospho-enrichment

Standard ProtiFi conditions are robust and different samples and experimental variations, especially targeted assays, may benefit from optimization.

Keep in mind:

- If maximal enzyme/substrate performance is desired, for example to obtain maximum liberation of a target MRM peptide, hold the sample type and amount constant, and test different enzyme:substrate ratios (e.g. 1:10, 1:50, and 1:100 at different times (e.g. 1- 3 hours).
- A typical protease:protein substrate ratio is 1:20 (w/w). A range of 1:10-1:100 can be implemented.
- A typical incubation period at 37 °C is 4- 16 hours (overnight). A temperature of 47 °C for a shorter duration (1- 2 hours) can also be used, enzyme dependent.
- Specifics for your process must be tested accordingly.
- Proteases:
 - Asp-N (Reference: Tarentino, A. L., et al. Molecular Cloning and Sequence Analysis of Flavastacin: An O-Glycosylated Prokaryotic Zinc Metalloendopeptidase. Archives of Biochemistry and Biophysics, 319 (1): 281-285 (1995). <https://doi.org/10.1006/abbi.1995.1293>)
 - Cleaves to the amino side of aspartic acid residues
 - Asp-N is a metalloprotease and requires small amounts of zinc to enhance activity
 - Maximal activity at pH 7.0- 8.0
 - Chymotrypsin (Reference: Kostka, V. and Carpenter, F. H. Inhibition of Chymotrypsin Activity in Crystalline Trypsin Preparations. JBC, 239 (6): 1799-1803 (1963). [https://doi.org/10.1016/S0021-9258\(18\)91261-5](https://doi.org/10.1016/S0021-9258(18)91261-5))
 - Cleaves to the carboxyl side of tyrosine, phenylalanine, tryptophan, and leucine residues
 - Treatment with N-alpha-p-tosyl-L-lysine chloromethyl ketone (TLCK) eliminates residual trypsin activity
 - Glu-C (Reference: Liu. S., et al. Mildly Acidic Conditions Eliminate Deamidation Artifact during Proteolysis: Digestion with Endoprotease Glu-C at pH 4.5. Amino Acids, 48 (4): 1059-1067 (2016). <https://doi.org/10.1007/s00726-015-2166-z>)
 - Cleaves the carboxyl side of glutamic acid
 - Activity and cleavage specificity is affected by buffer conditions